








Inventory Control Systems for vials

Model	Rack Arrangement	Storage phase	Number of triangular aluminum racks		Shelves per rack	Number of plastic dividers		Total number of 2ml vials
3K		Liquid	6		9	54		3 726
		Vapor	6		6	36		2 484
Model	Rack Arrangement	Storage phase	Number of aluminum racks		Shelves per rack	Number of cryoboxes		Total number of 2ml vials
			5x5	10x10		5x5	10x10	
10K		Liquid	4	7	13	52	91	10 400
		Vapor	4	7	11	44	77	8 800
24K		Liquid	6	17	13	78	221	24 050
		Vapor	6	17	11	66	187	20 350
38K		Liquid	6	28	13	78	364	38 350
		Vapor	6	28	11	66	308	32 450
LABS 20K		Vapor / Liquid	4	14	13	52	182	19 500
LABS 40K		Vapor / Liquid	8	30	13	104	390	41 600
LABS 80K		Vapor / Liquid	12	58	13	156	754	79 300

Racks: 10x 10 for 100vial boxes, 5x 5 for 25vial boxes

Other formats are available on demand