

XT Series



XT Series

The Extended Time Series (XT) Series of cryogenic refrigerators are designed for long-term storage of a variety of materials at cryogenic temperatures. This series also offers low profile models (XTL), with 5" canisters.



XT Series

| Models | XTL3 | XTL8 | XT10 | XT20 | XT34 | |
|--------------------------------------------|--------|-----------|--------|-----------|-----------|--------|
| Static Holding Time days (1) | 27 | 80 | 100 | 230 | 340 | |
| Working Time days (2) | 17 | 50 | 62 | 140 | 212 | |
| Evaporation Rate (1) liters/days | 0.11 | 0.10 | 0.10 | 0.09 | 0.10 | |
| Liquid Nitrogen Capacity - liters | 3 | 8 | 10 | 20.7 | 34 | |
| Weight Empty | - lbs. | 7.2 | 19.6 | 16.5 | 26 | 34.75 |
| | - kg. | 3.3 | 8.9 | 7.5 | 11.8 | 15.8 |
| Weight Full (3) | - lbs. | 12.5 | 33.8 | 34.3 | 62.8 | 95.3 |
| | - kg. | 5.7 | 15.4 | 15.6 | 28.6 | 43.3 |
| Neck Diameter | - in. | 2.0 | 2.0 | 2.0 | 2.0 | 2.0 |
| | - mm. | 51 | 51 | 51 | 51 | 51 |
| Overall Height | - in. | 17.2 | 19 | 23.8 | 25.8 | 26.3 |
| | - mm. | 437 | 483 | 597 | 655 | 668 |
| Overall Diameter | - in. | 7.6 | 15.6 | 11.4 | 15.6 | 18.8 |
| | - mm. | 193 | 396 | 290 | 396 | 478 |
| Number of Canisters | 6 | 6 | 6 | 6 | 6 | |
| Canister Dimensions (4) | - in. | 1.5x5 | 1.5x5 | 1.5x11 | 1.5x11 | 1.5x11 |
| | - mm. | 38x127 | 38x127 | 38x279 | 38x279 | 38x279 |
| Number of 1.2 ml & 2.0 ml vials (5/cane) | NA | NA | 150 | 150 | 150 | |
| Number of 1.2 ml & 2.0 ml vials (6/cane) | NA | NA | 180 | 180 | 180 | |
| Number of 1/2 cc straws (10/cane) | NA | NA | 540 | 540 | 540 | |
| Number of 1/2 cc straws - Bulk (1 level) | 750 | 750 | N/A | N/A | N/A | |
| Number of 1/2 cc straws - Bulk (2 levels) | NA | NA | 1500 | 1500 | 1500 | |
| Roller Base | N/A | R018-8C00 | N/A | R018-8C00 | R033-8C00 | |
| Low-Level Alarm | N/A | N/A | N/A | R033-8C15 | R033-8C15 | |

(1) Evaporation rate and static holding time are nominal. Actual rate may be affected by the nature of the contents, atmospheric conditions, container history, and manufacturing tolerances.

(2) Work time is an arbitrary, reference-only value to estimate container performance under the actual operating conditions.

(3) Without canisters.

(4) Canisters also available in 5.0 in. (127 mm) height for 10XT, 20XT and 34XT.

Features

- (1) New Rugged Construction – ribbed high strength aluminum body, magniformed necktube design, and more durable paint
- (2) Versatile storage system – convenient canister numbered index location ring and tapered internal spider
- (3) Superior vacuum performance with super insulation provides maximum holding times
- (4) Security – Accessory low-level alarm is available with remote monitor capability
- (5) Lockable lids
- (6) Mobility – Roller bases are available for some models



Taylor-Wharton

A legacy of excellence. A vision of promise.

www.taylorwharton.com
salesandservice@taylorwharton.com

The Americas

Taylor-Wharton Cryogenics LLC

4075 Hamilton Boulevard
Theodore, AL 36582
Phone: +251-443-8680
Phone: +800-898-2657

Australia

Taylor-Wharton Australia Pty. Ltd.

Unit 1/882 Leslie Drive
Albury, NSW 2640
Australia
Phone: +61 (0)2 6040-2533

Europe

Taylor-Wharton Germany GmbH.

Postfach 1470
D-25804 Husum
Phone: +49 (0)4841-9850

Taylor-Wharton Slovakia s.r.o.

Vstupný areál U.S. Steel
Q44 54 Košice
Slovakia
Phone: +421 (0)55-727-7157

Asia

**Taylor-Wharton (Beijing)
Cryogenic Equipment Co. Ltd.**

East End of Ping An Street
Xianghe Economic Development Zone
Xianghe City, Lang Fang Prefecture
Hebei Province, PRC. 065402
Phone: +86 (0)316-8219770

Taylor-Wharton Malaysia Sdn. Bhd.

Lot Nos. PT5073, 5076 & 5077
Jalan Jangur 28/43
Hicom Industrial Estate
40400 Shah Alam
Selangor Darul Ehsan
Malaysia
Phone: +60 (0)3-51913003

© 2012 Taylor-Wharton International. All rights reserved.